



## Quick Reference Guide

### Prevention & Management of Skin Injury from Personal Protective Equipment (PPE) – Addendum

Please note: Additional N95 Respirators used in Southern Health-Santé Sud have been tested for dressing compatibility – see table below.

<b>N95 Respirator Type</b>	<b>Possible Solutions</b>
Kimberly Clark N95 Regular	DO NOT use dressing
Kimberly Clark N95 Small	DO NOT use dressing
North 7130 Small	Thin non-bordered silicone foam
North 7130 Regular	DO NOT use dressing